Happy Intim® Rebirth Gel

LANGUAGE:

ENG.....	2
ESP.....	3
BUL.....	4
HRV.....	5
CES.....	6
DAN.....	7
NLD.....	8
EST.....	9
DEU.....	10
GRE.....	11
HUN.....	12
ITA.....	13
LAV.....	14
LIT.....	15
POL.....	16
POR.....	17
RON.....	18
SLO.....	19
SLV.....	20
SWE.....	21

ENG

DESCRIPTION

Non-cross-linked Hyaluronic Acid and protective buffer, class III medical device, made in conformity with Directive 93/42/EEC.

Hyaluronic Acid is from biotechnological, non animal origin. A medical device can only be injected by legally qualified physician in a legally authorized clinic.

PRESENTATION

Syringes containing 4 mL (0,1 fl Oz) / 5 mL (0,175 fl Oz) .

Boxes of 1, 2, 4 syringes.

COMPOSITIONS

Sodium Hyaluronate 60 mg/4mL / 75 mg/5mL
Protective buffer 16,8 mg/4mL / 21 mg/5 mL

INDICATIONS

Dermal implant of Hyaluronic acid for treatment of skin photo ageing and its consequences; treatment of dermal atrophy; HA is solubilized in a protective buffer solution composed by amino acids with hydration and antioxidant properties that stabilize the HA.

HOW TO USE

The medical device **HAPPY INTIM® REBIRTH GEL** is injected with or without papules in the dermis, using a 30G or a 32 G needle, with needle bevelled edge downwards. Nevertheless, the physician will use the method which best suits, and among others may be with: multiple longitudinal retro-tracing injections, fanning outwards, criss-cross technique, fish-bone technique or even transdermal techniques. Automated techniques for injections like mesotherapeutic guns may also be used. Any injection may leave small ecchymoses points or even haematomas, mainly on patients taking acetylsalicylic acid, anti-inflammatory or taking medication delaying coagulation. It does not often appear when the product is injected at a depth of more than 1 mm. Sensitive skins may benefit from the application of an anaesthetic cream prior to the treatment; said cream may cause redness or skin hypersensitivity. A too deep injection into the subcutaneous tissue does not produce the desired increase in tissue volume. The medical device **HAPPY INTIM® REBIRTH GEL** can eventually be mixed with other injectable solutions, depending on the needs of the physician. In this case it is the physician's responsibility to verify the compatibility and the suitability of the products he uses. Before implantation, physician should inform the patient about the possibilities of side effects (pain, redness, ecchymoses, stinging sensations and swelling, local inflammation, usually disappearing in 24 hours) and check that he is healthy. The medical device **HAPPY INTIM® REBIRTH GEL** must be used under appropriate aseptic conditions in an authorized clinic on healthy disinfected non-inflamed skin. Non aseptic use could induce infections. After the treatment, stick the traceability label on the patient's medical record.

The areas usually treated are:

- Internal face of arms
- Internal face of thighs
- Peri-umbilical area
- Abdomen
- Labia majora
- Perineum

RECOMMENDATIONS FOR PATIENTS

After the treatment and during 12 hours: avoid extreme temperatures, saunas and Hammam, direct exposure to sun or UV. Avoid apply on skin products other than those advised by the doctor. The physician must inform the patient of possible local reactions related to the implantation of this resorbable device: redness and/or local inflammation which disappears within 12-24 hours. Slight oedema or small ecchymoses may occur, which usually disappear within 24 or 48 hours. Should any inflammatory reaction or any other side effect has not disappeared within a week, the patient must inform the clinician about it immediately. The doctor will treat the patient via appropriate treatment adapted to the problem.

The doctor is requested to report any notable side effect related to the medical device **HAPPY INTIM® REBIRTH GEL** to **Skin Tech Pharma Group SL**, through the website www.skintechpharmagroup.com, by e-mail at info@skintechpharmagroup.com or by post.

CONTRAINDICATIONS

Allergy to any of the ingredients. Patients presenting any skin alteration, skin disease, infections or sequelae of streptococcal infections. Patients taking immunosuppressants, undergoing cortical therapy, with autoimmune disease history, patients with uncompensated diabetes, acute joint rheumatics, repetitive angina, and endocarditis. No studies are available for use during

pregnancy and breastfeeding or in case of treatment on children or minors under 18.

PRECAUTIONS

In case of allergic risk, a test should be done before the first injection. The test results must be attached to the patient's medical record. Do not re-use a syringe once opened. Do not use if the external protecting package has been seriously damaged or if you see some liquid inside of the box. Do not re-use or re-sterilize the unused leftover product; throw the leftover product away; throw the needle and the syringe into specific containers. Inform the patients under anticoagulant treatment about the risk of ecchymoses or haematomas due to injections. The patient should avoid taking acetylsalicylic acid or anti-inflammatory drug during the week preceding injections. Do not inject into the blood vessels, bones, tendons, muscles, ligaments or any organ different of skin. Prior to injecting, check the expiry date of syringes and needles; also check the integrity of the syringe seal. Do not use a product which is expired or looks damaged.

Caution: the outside of the presentation is not sterile. Once the presentation has been opened, use immediately. Do not heat, do not frost the presentation.


Inform the patient of product characteristics and possible side effects, precautions and contraindications.

There is no study available regarding cross tolerance or possible reaction between the medical device **HAPPY INTIM® REBIRTH GEL** and definitive filler products, even if it is not injected in the same plane or area. Should definitive implants exist, the medical device **HAPPY INTIM® REBIRTH GEL** may be used transdermally (protocol with Microneedling or electroporation, without injecting). Hyaluronic Acid precipitates when coming into contact with quaternary ammonium cations (such as Benzalkonium Chloride). Therefore, we do not recommend the medical device **HAPPY INTIM® REBIRTH GEL** coming into contact with said product.

STORAGE

Keep away from the light, under the temperature between 2 and 35 °C.

Do not freeze, do not heat.

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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Držite podalje od svjetla, na temperaturi između 2 i 35 ° C.

Ne zamrzavajte, ne zagrijavajte.

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2°C - 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restitilize
	Sterilized using steam or dry heat
	Temperature limitation 2°C - 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not sterilize
	Sterilized using steam or dry heat
	Temperature limitation 2°C - 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restitilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2°C — 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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
ARMAZENAMENTO

Manter afastado da luz, sob a temperatura entre 2 e 35 °C.

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restilize
	Sterilized using steam or dry heat
	Temperature limitation 2°C - 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirilizie
	Sterilized using steam or dry heat
	Temperature limitation 2°C - 35°C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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
Hranite stran od svetlobe, pri temperaturi med 2 in 35 °C.

Ne zamrzujte in ne segrevajte.

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Use by
	Do not use if package is damaged
	Do not reuse
	Do not restirize
	Sterilized using steam or dry heat
	Temperature limitation
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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