RRS[®] HA EYES

LANGUAGE:

ENG.....	2
SPA.....	3
HRV.....	4
CES.....	5
DAN.....	6
NLD.....	7
EST.....	8
HUN.....	9
LAV.....	10
LIT.....	11
POL.....	12
RON.....	13
SLO.....	14
SLV.....	15
SWE.....	16

ENG

DESCRIPTION

RRS® HA EYES is a sterile medical device, class III, made in conformity with Directive 93/42/EEC, intradermal implant that contains a non-cross-linked, resorbable hyaluronic acid associated to protective buffer that helps to prolong the effect of the dermal implant, specially formulated for ageing skin around the eyes. Hyaluronic acid is from biotechnological, non-animal origin. A medical device can be injected only by a legally authorised physician in a legally authorised clinic.

PRESENTATION

Ampoules containing 1.5 ml (0,05 fl Oz.); boxes of 2 amp, 12 amp, 24 amp.

COMPOSITION

Sodium Hyaluronate 7,5 mg / 1,5 ml

Protective buffer 40,95 mg / 1,5 ml

INDICATIONS

Dermal implant for treatment of skin photo ageing and its consequences or skin atrophy especially of the area around the eyes.

HOW TO USE

RRS® HA EYES is injected in the superficial dermis, using a 32 G needle. We recommend the use of a simple multipuncture, very superficial protocol, including the temporal area and injecting a very little quantity at each point in order to limit frequent bruising and possible stinging sensations on this very sensitive area. A maximum of 0,75 ml is sufficient for each side. Automated techniques for injections like mesotherapy guns may also be used. Any injection may leave small ecchymoses points or even haematomas, mainly on patients taking acetylsalicylic acid, anti-inflammatory or taking medication delaying coagulation. Sensitive skins may benefit from the application of an anaesthetic cream prior to the treatment; said cream may cause redness or skin hypersensitivity.

RRS® HA EYES can eventually be mixed with other injectable solutions, depending on the needs of the physician. In this case it is the physician's responsibility to verify the compatibility and the suitability of the products he uses. Before implantation the physician should inform the patient about the possibilities of side effects (pain, redness, ecchymoses, stinging sensations, swelling, local inflammation, usually disappearing in 24 hours) and check that he is healthy. RRS® HA EYES must be used under appropriate aseptic conditions in an authorised clinic on healthy, disinfected, non-inflamed skin. Non aseptic use could induce infections. After the treatment, stick the traceability label on the patient's medical record.

The areas usually treated are:

- Area around the eyes, up to temporal area
- Bags under the eyes, circles, dark circles

RECOMMENDATIONS FOR PATIENTS

After the treatment and during 12 hours: avoid extreme temperatures, saunas and Hammam, direct exposure to sun or UV. Avoid to make-up or apply on skin products other than those advised by the doctor. The physician must inform the patient of possible local reactions related to the implantation of this resorbable device: redness and/or local inflammation which disappears within 12-24 hours. Slight oedema or small ecchymoses may occur, which usually disappear within 24 or 48 hours. Should any inflammatory reaction or any other side effect has not disappeared within a week, the patient must inform the clinician about it immediately. The doctor will treat the patient via appropriate treatment adapted to the problem.

The doctor is requested to report any notable side effect related to the RRS® HA EYES to Skin Tech Pharma Group SL. through the website www.skintechpharmagroup.com.

com, by e-mail at info@skintechpharmagroup.com or by post.

CONTRAINDICATIONS

Allergy to any of the ingredients. Patients presenting any skin alteration, skin disease, infections, or sequelae of streptococcal infections. Patients taking immunosuppressants, undergoing cortical therapy, with autoimmune disease history, patients with uncompensated diabetes, acute joint rheumatics, repetitive angina, and endocarditis. No studies are available for use during pregnancy and breastfeeding or in case of treatment on children or minors under 18.

PRECAUTIONS

In case of allergic risk, a test should be done before the first injection. The test results must be attached to the patient's medical record. Do not re-use ampoule once opened. Do not use if the external protecting package has been seriously damaged or if you see some liquid inside of the box. Do not re-use or re-sterilise the unused leftover product; throw the leftover product away; throw the needle and the ampoule into specific containers.

Inform the patients under anticoagulant treatment about the risk of ecchymoses or haematomas due to injections. The patient should avoid taking acetylsalicylic acid or anti-inflammatory drug during the week preceding injections. Avoid direct contact with eyes. Do not inject into the blood vessels,

bones, tendons, muscles, ligaments or any organ different of skin. Prior to injecting, check the expiry date of ampoules, syringes and needles; also check the integrity of the ampoule seal. Do not use a product which is expired or looks damaged.

Caution: the outside of the presentation is not sterile. Once the presentation has been opened, use immediately. Do not heat, do not frost the presentation. Inform the patient of product characteristics and possible side effects, precautions, and contraindications.

There is no study available regarding cross tolerance or possible reaction between RRS® HA EYES and definitive filler products, even if it is not injected in the same plane or area. Should definitive implants exist, RRS® HA EYES may be used transdermally (protocol with AD ROLL TD® or ExcellDerm® Pro, without injecting). Hyaluronic acid precipitates when coming into contact with quaternary ammonium cations (such as benzalkonium chloride), therefore we do not recommend RRS® HA EYES coming into contact with said product.

STORAGE

Keep away from the light, under the temperature between 2 and 35 oC.

Do not freeze, do not heat..

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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No congelar, no calentar.

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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Držite podalje od svjetla, na temperaturi između 2 i 35 ° C.

Ne zamrzavajte, ne zagrijavajte.

	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C — 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C – 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: SKIN TECH PHARMA GROUP, S.L. C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilizze
	Sterilized using steam or dry heat
	Temperature limitation
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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