# RRS<sup>®</sup> HA STRIMATRIX

## LANGUAGE:

ENG.....	2
SPA.....	3
HRV.....	5
CES.....	7
DAN.....	8
NLD.....	9
EST.....	10
HUN.....	11
LAV.....	12
LIT.....	13
POL.....	14
RON.....	16
SLO.....	17
SLV.....	18
SWE.....	19

## ENG

### DESCRIPTION

**RRS® HA STRIMATRIX** is a sterile medical device, class III, made in conformity with Directive 93/42/EEC, intradermal implant that contains a non-cross-linked, resorbable hyaluronic acid associated to protective buffer that helps to prolong the effect of the dermal implant, specifically formulated for skins with striae distensae and flaccidity. Hyaluronic acid is from biotechnological non animal origin. A medical device only can be injected by a legally authorised physician in a legally authorised clinic.

### PRESENTATION

Vials containing 10 ml (0,35 Fl.oz.); boxes of 2 vials, 6 vials, 12 vials.

### COMPOSITION

Sodium Hyaluronate	20 mg / 10 ml
Protective buffer	659,1 mg / 10 ml

### INDICATIONS

**RRS® HA STRIMATRIX** Dermal implant of non cross-linked hyaluronic acid for the treatment of stretch marks, to fill the depressions caused by linear skin ruptures of stretch marks. The activity of the implant is sustained by the associated protective buffer which allows to regenerate the dermal matrix, atrophic depressions, and flaccidity of the dermis at the level of stretch marks.

### HOW TO USE

**RRS® HA STRIMATRIX** is injected as superficially as possible in the dermis, with a 32 G needle, using the technique of multiple superficial longitudinal retro-tracing injections. The physician nevertheless can use the method which best suits, and among others may be with: multiple punctures, fanning outwards, fish bone technique. Automated techniques for injections like mesotherapy guns may also be used.

Any injection may leave small ecchymoses points or even haematomas, mainly on patients taking acetyl salicylic acid, anti inflammatory in patients taking medication delaying coagulation. Sensitive skins may benefit from application of an anaesthetic cream prior to the treatment; said cream may cause redness or skin hypersensitivity. **RRS® HA STRIMATRIX** can eventually be mixed with other injectable solutions, depending on the needs of the physician. In this case it is the physician's responsibility. to verify the compatibility and the suitability of the products he uses.

Before implantation, the physician should inform the patient about the possibilities of side effects (pain, redness, ecchymoses, stinging sensations, swelling, local inflammation, usually disappearing in 24 hours) and check that he is healthy. **RRS® HA STRIMATRIX** must be used under appropriate aseptic conditions in an authorised clinic on healthy, disinfected, non-inflamed skin. Non aseptic use could induce infections. After the treatment, stick the traceability label on the patient's medical record.

The areas usually treated are:

- Abdomen
- Buttocks
- Internal face of arms
- Internal face of thighs
- Arms, shoulders

### RECOMMENDATIONS FOR PATIENTS

After the treatment and during 12 hours: avoid extreme temperatures, saunas and Hammam, direct exposure to sun or UV. Avoid to apply on skin products other than those advised by the doctor. The physician must inform the patient of possible local reactions related to the implantation of this resorbable device: redness and/ or local inflammation which disappears within 12-24 hours. Slight oedema or small ecchymoses may occur, which usually disappear within 24 or 48 hours. Should any inflammatory reaction or any other side effect has not disappeared within a week, the patient must inform the clinician about it immediately. The doctor will treat the patient via appropriate treatment adapted to the problem. The doctor is requested to report any notable side effect related to the **RRS® HA STRIMATRIX** to **Skin Tech Pharma Group SL** through the website [www.skintechpharmagroup.com](http://www.skintechpharmagroup.com), by e-mail at [info@skintechpharmagroup.com](mailto:info@skintechpharmagroup.com) or by post.

### CONTRAINDICATIONS

Allergy to any of the ingredients. Patients presenting any skin alteration, skin

disease, infections or sequelae of streptococcal infections. Patients taking immunosuppressants, undergoing cortical therapy, with autoimmune disease history, patients with uncompensated diabetes, acute joint rheumatics, repetitive angina, and endocarditis. No studies are available for use during pregnancy and breastfeeding or in case of treatment on children or minors under 18.

### PRECAUTIONS


In case of allergic risk, a test should be done before the first injection. The test results must be attached to the patient's medical record. Do not re-use vial once opened. Do not use if the external protecting package has been seriously damaged or if you see some liquid inside of the box. Do not re-use or re-sterilise the unused leftover product; throw the leftover product away; throw the needle and the vial into specific containers. Inform the patients under anticoagulant treatment about the risk of ecchymoses or haematomas due to injections. The patient should avoid taking acetylsalicylic acid or anti inflammatory drug during the week preceding injections. Do not inject into blood vessels, bones, tendons, muscles, ligaments or any organ different of skin. Prior to injecting, check the expiry date of vials, syringes and needles; also check the integrity of the vial seal. Do not use a product which is expired or looks damaged.

Caution: the outside of the presentation is not sterile. Once the presentation has been opened, use immediately. Do not heat, do not frost the presentation. Inform the patient of product characteristics and possible side effects, precautions and contraindications.

There is no study available regarding cross tolerance or possible reaction between **RRS® HA STRIMATRIX** and definitive filler products, even if it is not injected in the same plane or area. Should definitive implants exist, **RRS® HA STRIMATRIX** may be used transdermally (protocol with **AD ROLL TD®** or **ExcellDerm® Pro**, without injecting.) Hyaluronic acid precipitates when coming into contact with quaternary ammonium cations (such as benzalkonium chloride), therefore we do not recommend **RRS® HA STRIMATRIX** coming into contact with said product.

### STORAGE

Keep away from the light, under the temperature between 2 and 35 °C. Do not freeze, do not heat.

	Manufacturer / Fabricant / Fabricante: <b>SKIN TECH PHARMA GROUP, S.L.</b> C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: <b>SKIN TECH PHARMA GROUP, S.L.</b> C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) <a href="http://skintechpharmagroup.com">skintechpharmagroup.com</a> <a href="mailto:info@skintechpharmagroup.com">info@skintechpharmagroup.com</a>
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: <b>SKIN TECH PHARMA GROUP, S.L.</b> C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) skintechpharmagroup.com info@skintechpharmagroup.com
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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	Manufacturer / Fabricant / Fabricante: <b>SKIN TECH PHARMA GROUP, S.L.</b> C/ Pla de l'Estany, 29 17486 Castelló d'Empúries · Girona (SPAIN) <a href="http://skintechpharmagroup.com">skintechpharmagroup.com</a> <a href="mailto:info@skintechpharmagroup.com">info@skintechpharmagroup.com</a>
	Date of manufacture
	Batch code
	Do not use if package is damaged
	Do not reuse
	Do not restirilize
	Sterilized using steam or dry heat
	Temperature limitation 2 °C - 35 °C
	Does not contain latex
	Keep away from sunlight
	Keep dry
	Notified body number

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